CE 0483

Orthodontic Bonding System

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Dear Customer

Thank you for having chosen a quality product from Dentaureum.

In order to use this product at its best for you and your patients, it is important to study and follow these directions for use carefully.

The written instructions cannot cover all eventualities during operation. For this reason our Hotline is available to answer any other questions and ideas that may arise.

Due to a constant development we recommend, even when you use the same products frequently, to study the enclosed latest issue of the instructions for use or refer to our website at www.dentaureum.de.

Manufacturer

Dentaureum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

Contents

1. Safety references – please read carefully!	20
2. General instructions for adhesive technique – please read carefully!	21
3. Tips for debonding – please read carefully!	24
4. NO-MIX one-component bracket adhesive	26
5. MIX two-component bracket adhesive	29
6. With reference to our quality standards	32
7. Explanation of symbols used on the label.....	33

1. Safety references – please read carefully!

1.1. NO-MIX and MIX bracket adhesives:

- Medical products, store in a safe place, recommended storage temperature 5–15 °C / 41–59 °F in the refrigerator. Process at room temperature, remove from refrigerator at least 30 minutes before work begins. This assures that the adhesive reaches room temperature (max. 25 °C/77 °F).
- Close immediately after use. Do not use after expiration date has expired.
- Keep away from children.

1.2. Etching Gel

Etching gel contains 37.4 % O-phosphoric acid. Avoid contact with dentin, soft tissue, eyes and skin. If contact occurs on these surfaces, flush clean with water. In case of contact with eyes, get medical attention.

1.3. Pastes and sealant resins, MIX System

Avoid prolonged contact with soft tissue and skin, especially for people who are allergic to acrylics. Both components, the sealant resin, and the pastes, should not come in contact with one another to avoid contamination.

Patients, who are allergic to Acrylics or Dibenzoylperoxide must not be treated with these adhesives!

1.4. Activator and adhesive paste, NO-MIX-System

- Avoid prolonged and repeated contact with skin. In case of contact, wash skin thoroughly with soap and water. Uncured materials may cause contact dermatitis. Apply in a well-ventilated area.
- Avoid contamination of the paste with the activator.

Patients, who are allergic to Acrylics or Dibenzoylperoxide must not be treated with these adhesives!

1.5. Primer for plastic brackets

Avoid prolonged or repeated contact with the skin. If it occurs, wash skin thoroughly with soap and water. Contact may cause skin irritation. Apply in a well-ventilated area.

Caution: Do not expose primer to sunlight. Do not store over 25 °C / 77 °F.

1.6. Ceramic brackets

When bonding ceramic brackets, it is essential to avoid contact with the antagonist due to the risk of abrasion. If necessary, a palatal bite raise can be used.

2. General instructions for adhesive technique – please read carefully!

2.1. Checking the tooth surface

Check each tooth carefully for any existing enamel damage (cracks, lesions or similar). Such predamages may lead to enamel fractures when debonding brackets. The patient must be informed accordingly prior to the treatment.

2.2. Cleaning the teeth

The tooth surfaces to be bonded must first be thoroughly cleaned. Small brushes or rubber cups with special prophylactic paste or a pumice and water mix can be used. It is essential that the cleaning agent contains no fluoride, oil or grease since these would interfere with future etching procedures. Slow handpiece speed should be used.

2.3. Dry field

Successful adhesive procedures require a dry, clean and visible working field. With most patients, an acceptable dry field is achieved with cotton rolls. If extensive bonding procedures are to be undertaken, the use of a rubber dam is advisable in order to achieve perfect dryness. This is particularly true when a patient has a heavy saliva flow or other unfavorable working conditions exist.

2.4. Enamel etching

The enamel surface is etched with the etching gel. A characteristic surface relief with 20 µm deep crevices is obtained through selective dissolution of the enamel prism. This results in an increased surface area due to roughening and permits intimate mechanical interlocking with the adhesive. Use only a dabbing action, do not rub the etching gel onto the enamel.

The required etching time is at least 15 seconds. In special cases e.g. in cases of heavily fluoridated teeth, the etching time must be extended by at least 10 seconds.

Avoid etching gel contact with the gingiva, skin, eyes and dentin.

2.5. Rinsing off etching gel residue and drying the etched surfaces

Etching gel residues must be rinsed off with a strong jet of water (at least 5 seconds for each tooth). If a combined air and water syringe is used to dry the teeth, it is essential that air from the compressor contain no oil. In case of doubt, special ether-based drying agents or a modified dryer should be used to dry the teeth. The clinical criterion of satisfactory etching is a uniform frosty white appearance of etched teeth.

The etched surfaces should not be touched and must be kept dry and clean!

2.6. Selection, bonding and adaptation of brackets

Metal, ceramic or plastic brackets may be used. The curvature of the base must correspond to that of the tooth surface to be bonded.

2.6.1. Bonding with plastic brackets

When using the NO-MIX or MIX-adhesive, the special **primer (REF 163-500-00)** for Elegance® plastic brackets must be applied. The primer supplies a chemical bond between the plastic bracket and the adhesive.

Procedure: with a small brush, located in the lid of the jar, apply a thin and even layer of primer to the base of the bracket. For best bonding results we recommend applying the adhesive within 30 seconds. The maximum drying time of 60 seconds should not be exceeded. Continue following bonding information from page 26 onwards.

Attention: Do not expose primer to direct sunlight.

Maximum storage temperature 25 °C / 77 °F.

These specifications pertain only to Elegance® plastic brackets. When using the Orthodontic Bonding System for plastic brackets from other manufacturers, please follow their instructions.

2.6.2. Bonding with ceramic brackets

Thanks to the laser-structured base, discovery® pearl ceramic brackets can be used like metal brackets. The laser-structure secures the mechanical bond to the NO-MIX and MIX adhesives.

However, the base of the discovery® pearl brackets should not be touched with your fingers or contaminated in any way, as this could otherwise detrimentally affect the bond. The silanization process on the base of our Fascination®, Fascination® 2 and Jewels ceramic brackets in connection with the Orthodontic Bonding System supplies a very good chemical bond. We recommend the MIX-adhesive when using Fascination® 2 ceramic brackets.

When using the Orthodontic Bonding System for ceramic brackets from other manufacturers, please follow their instructions.

2.6.3. Bonding to crowns and fillings

Bonding to ceramic crowns requires the use of a ceramic primer such as Rely X™ Ceramic Primer from 3M™ Espe™.

Bonding to **acrylic crowns** does not require any special conditioning; if necessary, the special primer (REF 163-500-00) for Elegance® plastic brackets can be applied to the acrylic crown.

Bonding to **metal crowns** requires roughening, e.g. using a MicroEtcher. Then apply, e.g. Assure Bonding Resin from Reliance Orthodontic Products or Rely X™ Ceramic Primer from 3M™ Espe™.

Do not etch the crowns!

Adhere to the instructions for use of the materials used!

2.7. Bracket positioning and ligation of the arch wire

Apply adhesive to the bracket and use a clamping tweezer (e.g. REF 025-277-00), seat the bracket on the prepared enamel surface exerting light pressure. Within 15 seconds, recheck the position of the bracket and adjust if necessary. Press the bracket onto the tooth surface so that any excess adhesive squeezes out from the edges. If the correct amount of adhesive has been applied to the bracket, there will be little excess to be removed. Any excess can be removed, also in the interdental areas, with a scaler (e.g. REF 027-349-00) or an instrument for the bonding technique (e.g. REF 023-275-00). Once polymerization has started the bracket should not be touched and the adhesive should set without being disturbed. Arch wire ligation can begin 10 minutes after last bracket has been seated.

3. Tips for debonding – please read carefully!

Debonding has to be performed with great care!

Adhesive residue can be removed with a scaler (e.g. REF 027-349-00) or with a carbide bur (REF 123-603-00 or 123-604-00) at low speed (recommended speed 40000 min⁻¹), in a dry working area and with a continuous air stream. Then the enamel should be polished, using a rubber cup and a pumice-water mixture.

A fluoride treatment is advisable after final polishing!

Alternatively adhesive residue removers (REF 123-700-00 to 123-703-00) can be used (recommended speed 5000 – max. 10000 min⁻¹ using spray cooling).

3.1. Metal and plastic brackets

We recommend using Weingart pliers (REF 003-120-00) for debonding our metal and plastic brackets.

Grip the metal brackets mesially and distally or diagonally on the wings. Lightly press the jaws of the pliers together and use the predetermined bending point of the discovery®/discovery® smart brackets and the curvature of the equilibrium® brackets to carefully loosen the base from the tooth.

To debond Elegance® plastic brackets grip the brackets mesially and distally on the base. Then remove the brackets from the tooth **without pulling** by twisting them briefly to the left or right.

Please adhere to the debonding instructions in the respective instructions for use of the brackets.

3.2. **discovery® pearl, Fascination®, Fascination® 2 and Jewels ceramic brackets discovery® pearl**

Remove arch wire from the bracket. For debonding the discovery® pearl brackets we recommend the special discovery® pearl debonding instrument, REF 019-001-00, which grips the body of the bracket directly at the base. The bracket is gently removed from the bracket adhesive by slightly tipping it sideways in the direction of the slot and archwire.

Alternatively, the straight bracket removal pliers (REF 004-346-00) can be used for debonding the discovery® pearl brackets. Place the blades of the bracket removal pliers mesio-distally on the bracket base, directly on the surface of the tooth. The bracket is gently removed from the tooth by slightly twisting the pliers. The archwire can be removed beforehand – though it is not essential.

Despite their small, delicate shape there is no risk of fracture with discovery® pearl brackets provided they are used correctly.

Fascination®, Fascination® 2 and Jewels ceramic brackets

For the mechanical debonding we recommend the Weingart pliers (REF 003-120-00).

Remove arch wire from the bracket. To avoid damage to the bracket while debonding, place a blue separating ring (e.g. REF 774-200-01) around the bracket base. Grasp the bracket base mesially/distally with a Weingart plier (REF 003-120-00) at the bracket base directly at the tooth surface. Without any pressure, but with a slight twist, remove the bracket. **(DO NOT PULL).**

Hints:

Use of other instruments can result in bracket fracture. If the bracket fractures during **incorrect** debonding, we recommend carefully removing the fractured segments and the adhesive from the enamel using bracket debonding pliers (e.g. REF 003-349-00) or a scaler (e.g. REF 027-349-00).

Caution: Do not damage tooth enamel, danger of enamel fracture. Use proper evacuation procedures, protect everyone's eyes within the treatment area!

Attention: When using ceramic brackets from other manufacturers, please follow their instructions for bracket removal.

4. NO-MIX one-component bracket adhesive

4.1. Indications and characteristics

NO-MIX is a paste/liquid adhesive for the direct bonding of metal, ceramic and plastic brackets. The polymerization on the tooth is initiated on contact by pressing the bracket onto the tooth, so that time consuming mixing can be eliminated. When debonding the brackets separation is widely shifted to the interface between bracket and adhesive. This avoids damage to the enamel. The bracket adhesive in syringe form should be used sparingly and can be accurately applied. Excess adhesive around the base is usually unpolymerized and easy to remove.

4.2. Hints for use

Preparation of the teeth

- Cleaning
- Drying
- Etching
- Rinsing off etching residue and drying

Selection and adaptation of the brackets

For detailed information about the different applications refer to pages 20 – 26.

4.3. Activator application

- Use a disposable brush to apply a thin coating of the activator to both the etched and dried tooth surface and the bracket base. Excess activator will cause a dilution of the adhesive, resulting in bracket drifting and weakened bonds.

4.4. Application of the adhesive paste and ligation of the arch wire

- Apply adhesive paste with spatula or syringe directly onto the bracket base. With the syringe, the adhesive paste can be accurately dosed.
- Avoid contamination between the adhesive paste and activator.
- The bracket is positioned onto the tooth applying light pressure and can still be adjusted within approximately 15 seconds. After the final position is determined press the bracket firmly onto the tooth. The excess adhesive can be removed with a scaler (e.g. REF 027-349-00) or an instrument for the bonding technique (e.g. REF 023-275-00).
- **The arch wire can be placed 10 minutes after the last bracket has been bonded** (refer to page 24 (2.7)).

4.5. Ingredients

Adhesive paste:

Bis-GMA
Triethylenglycoldimethacrylate
Barium Glass
Pyrogenic Silica
Benzoyl Peroxide, stabilizer

Etching Liquid:

O-phosphoric acid

Activator:

Bis-GMA
Triethylenglycoldimethacrylate
Polymethylmethacrylate
Pyrogenic Silica
Amine starter, stabilizer

4.6. Storage and Shelf Life

- Please check expiration date on the label. When out of date, please do not use. Recommended storage temperature: 5 – 15 °C / 41 – 59 °F in the refrigerator.
- Remove package from refrigerator at least 30 minutes before use to guarantee that material reaches room temperature (max. 25 °C / 77 °F).
- Close immediately after use.

4.7. Disposal

Adhesive and adhesive residue can be disposed of with ordinary waste removal.

4.8. Availability

NO-MIX one-component bracket adhesive in jars

Contents (set)

1 x bracket adhesive 8 g, 1 x activator 10 ml,
1 x etching gel 10 ml, 1 x brush holder,
50 x disposable brushes, 1 x spatula,
1 x mixing pad

Refill packs / Accessories

A = bracket adhesive in jar
B = activator
C = etching gel in bottle
D = disposable brushes
E = spatula
F = brush holder
G = mixing pad
H = primer for plastic brackets
(not included in the set)

	163-406-00**	1 set
for approx. 800 brackets (40 treatment cases)		
A	163-413-00**	8 g jar
B	163-414-00**	10 ml
C	163-415-00**	10 ml
D	163-417-00*	50 pieces
E	163-418-00	50 pieces
F	163-419-00*	5 pieces
G	163-416-00	2 pieces
H	163-500-00**	10 ml

NO-MIX one-component bracket adhesive in syringes

Contents (set)

3 x bracket adhesive in syringes each 2.5 g,
1 x activator 10 ml, 2 x etching gel in syringes each 2.5 g,
6 x disposable tubes for etching gel syringes, 1 x brush holder, 50 x disposable brushes, 1 x spatula, 1 x mixing pad

Refill packs / Accessories

A = bracket adhesive in syringes
B = activator
C = etching gel in syringes
D = disposable brushes
E = spatula
F = brush holder
G = disposable plastic tubes for syringes
H = mixing pad
I = primer for plastic brackets
(not included in the set)

	163-407-00**	1 set
for approx. 600 brackets (30 treatment cases)		
A	163-421-00**	3 x 2.5 g syringes
B	163-414-00**	10 ml
C	163-412-00**	2 x 2.5 g syringes
D	163-417-00*	50 pieces
E	163-418-00	50 pieces
F	163-419-00*	5 pieces
G	163-422-00*	10 pieces
H	163-416-00	2 pieces
I	163-500-00**	10 ml

* CE / ** CE 0483

5. MIX two-component bracket adhesive

5.1. Indications and characteristics

Mix is a bracket adhesive consisting of an activating paste (A) and catalyst paste (B) which are mixed in a ratio of 1 : 1. The material is of a very soft consistency with sufficient flowing ability to achieve a good bonding strength, however, with enough viscosity to allow the brackets to be positioned without drifting.

5.2. Hints for use

Preparation of the teeth

- Cleaning
- Drying
- Etching
- Rinsing off etching residue and drying

Selection and adaptation of the brackets

For detailed information about the different applications refer to pages 20 – 26.

5.3. Application

5.3.1. Use of sealant resins A+B

Mix resins A + B (small red and white capped bottles) in a ratio of 1 : 1 for 5 – 10 seconds on mixing pad. Apply mixture in a thin coating, using application sponges on the etched areas. Suggested instruments: Mosquito forceps (REF 000-730-00) and/or dental tweezers. The sealant resin flows into the retention grooves and builds a mechanical bond with the etched surface.

5.3.2. Use of adhesives A+B

Before use, stir pastes A+B well. Mix pastes A+B (red and white capped jars) in a ratio of 1:1 with the mixing spatula on the mixing pad for approximately 20 seconds.

When removing the pastes from the jars, please be sure to use opposite ends of the mixing spatula in order to avoid contamination between pastes A+B. During mixing, make sure that no air bubbles are present as they will affect the cure of the paste.

5.3.3. Application of bonding agent on bracket base and ligation of the arch wire

The mixture of pastes is applied in a **thin** layer on the bracket base.

The bracket is positioned onto the tooth applying light pressure and can still be adjusted within approximately 15 seconds. After the final position is determined press the bracket firmly onto the tooth. The excess adhesive can be removed with a scaler (e.g. REF 027-349-00) or an instrument for the bonding technique (e.g. REF 023-275-00).

5.3.4. The arch wire may be placed 10 minutes after the last bracket has been bonded

(refer to page 24 (2.7))

5.4. Ingredients

Paste A:

Bis-GMA
Triethylenglycoldimethacrylate
Barium Glass
Quartz
Pyrogenic Silica
Amine starter, stabilizers, pigments

Paste B:

Bis-GMA
Triethylenglycoldimethacrylate
Barium Glass
Quartz
Pyrogenic Silica
Benzoyl peroxide, stabilizers

Etching gel:

O-phosphoric acid

Sealant resin A:

Bis-GMA
Triethylenglycoldimethacrylate
Amine starter
Stabilizers

Sealant resin B:

Bis-GMA
Triethylenglycoldimethacrylate
Benzoyl peroxide
Stabilizers

5.5. Storage and Shelf Life

- Please check expiration date on the label. When out of date, please do not use. Recommended storage temperature: 5–15 °C / 41–59 °F in the refrigerator.
- Remove package from refrigerator at least 30 minutes before use to guarantee that material reaches room temperature (max. 25 °C / 77 °F).
- Close immediately after use.

5.6. Disposal

Adhesive and adhesive residue can be disposed of with ordinary waste removal.

5.7. Availability

MIX two component bracket adhesive in jars

Contents (set)

1 x paste A 18 g, 1 x paste B 18 g,
1 x sealant resin A 3 ml, 1 x sealant
resin B 3 ml, 1 x etching gel 10 ml,
50 x spatulas, 2 x mixing pads,
100 x application sponges

Refill packs / accessories

A = paste A in jar
B = paste B in jar
C = sealant resin A in bottle
D = sealant resin B in bottle
E = etching gel in bottle
F = mixing pad
G = spatula
H = application sponges
I = primer for plastic brackets
(not included in the set)

163-405-00**	1 set
---------------------	-------

for approx. 1400 brackets
(70 treatment cases)

A	163-408-00**	18 g jar
B	163-409-00**	18 g jar
C	163-410-00**	3 ml
D	163-411-00**	3 ml
E	163-415-00**	10 ml
F	163-416-00	2 pieces
G	163-418-00*	50 pieces
H	163-420-00*	100 pieces
I	163-500-00**	10 ml

* CE / ** CE 0483

6. With reference to our quality standards

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

Please refer to the label.

**Additional information can be found at www.dentaurum.de
(Explanation of the symbols REF 989-313-00).**

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